Prior to use, carefully read the instructions for use

Cytrans™ Elashield

RESORBABLE P(LA/CL) (L-LACTIDE-ε-CAPROLACTONE COPOLYMER) BILAYER MEMBRANE

Cytrans™ Elashield is intended to be used as a barrier membrane in the regeneration of bone defects. Cytrans™ Elashield is used in conjunction with bone grafting procedures to prevent soft tissue infiltration into the grafted space and enhance bone formation.

INDICATIONS FOR USE

Cytrans™ Elashield is intended for the following uses as a space-making barrier:

- Guided bone regeneration in oral surgery
- Alveolar ridge augmentation
- Guided tissue regeneration in periodontal defects
- Extraction socket site preservation
- Sinus lifts
- Immediate implant placement at time of extraction or delayed placement
- Treatment of associated cystic defect
- The containment of bone grafting materials

CONTRAINDICATIONS

1. Patients receiving kidney dialysis
2. Patients with blood disorders
3. Patients suspected of having collagen disease
4. Patient receiving steroids
5. Diabetic patients who are either untreated or showing no improvement with drug treatment
6. Patients with fever or other infection
7. Immunodeficient patients
8. Patients with serious comorbidity (heart disease, etc.)
9. Pregnant women, lactating patients, and patients who may be pregnant
10. Patients receiving radiotherapy or having a history thereof
11. Patients being treated with bisphosphonates or having a history thereof
12. Patients with hypertension or hypotension
13. Asthma patients
14. Patients with rheumatic fever
15. Patients with hyperthyroidism
16. Patients with liver disease or thought to have an abnormal liver function
17. Patients judged by doctors or dentists as having inappropriate lifestyles and hygiene
18. Patients with a history of sensitivity to drugs, food, accessories, or chemical substances: the medical histories of these patients should be taken, and the product should be used with caution even in the absence of any history of sensitivity to the product or similar items.

PRECAUTIONS

I. Precautions for use

1. Do not re-use

2. Do not re-sterilize

II. Additional basic precautions

1. GBR surgery using the product should be performed by dentists familiar with this type of surgery.
2. To prevent contamination and infection when using the product, handle it with sterilized instruments. At the same time, be careful to ensure that the material does not come into contact with unsterilized objects such as teeth or saliva around affected areas.
3. Do not use the product if the package is damaged, as its sterility may be compromised.
4. To prevent creating perforations in the product avoid using tweezers, etc. with sharp tips.
5. Remove corners when trimming the product to prevent mucosal irritation.
6. When covering bone defect using the product, first completely remove any infected tissue or soft tissue from the defect areas.
7. Peel off the protective film before using the product.
8. Cover an area 2 mm or more from the margin of the bone defect with the product to prevent gingival tissue ingrowth from outside of the defect.
9. Avoid excessive stretching.
10. Note that if the glossy dense layers (solid layers) of the product are pressed together firmly, they may stick to one another and impair usability.
11. When suturing, the product is to be completely covered with a gingival flap, and care must be taken so that the product itself is not exposed or moved.
12. Do not use for purposes other than those stated in [Indication for use].
13. As this is a sterile product, it should be used immediately after opening the container and should be used only once.
14. Consider the use of a systemic antibiotic in those patients who may have reduced healing capacity.

DIRECTIONS FOR USE

1. Prior to use of the product, take measures to ensure hygiene such as oral hygiene and bacterial infection control, including appropriate instructions for patients.
2. Surgically expose bone defects via normal surgical procedures.
3. Fill the bone defect with the bone graft material.
4. Confirm the product size required, check the packaging material for damage and aseptically remove the product.
5. Use sterilized scissors to trim the product to a shape and size sufficient to cover the defect. After trimming, peel the product from the protective film. In addition, to facilitate maneuverability, the product can be wetted with the patient's blood or sterile saline before use.
6. Cover the bone defect. When doing so, the porous layer should face the bone defect, and the dense layer (solid layer) should face the gingival flap. You can identify each layer because the solid layer surface is obviously glossier in appearance than the porous layer.
7. The flap has to be repositioned with a tension-free closure. Make releasing incisions in the gingival flap as necessary.
8. Perform regular post-surgery examinations. When performing surgical procedures such as implants at the same location after a few weeks to a few months of follow-up, and at the discretion of a dentist after bone formation using this product, check, via inspection or palpation, for residue of the product and whether the bone surface is exposed. If any residue of the product is found, perform curettage removal and clean using periosteal elevator.

STORAGE AND HANDLING

I. Storage

Recommended for optimal performance, store in a cool and dark place (4-25℃ / 39-77°F).

II. Period of Use

This product is for single use only and should be disposed of after each single treatment of a single patient. If mixed with patient blood, the discarded material is considered biohazardous waste.

III. Expiration

The product should be used by the expiry date indicated on the packaging*.

*For example, " 2024-03-13" means the expiry date is March 13, 2024.

PACKAGE

Sizes

15 mm x 25 mm

25 mm x 25 mm

30 mm x 40 mm

Caution: U.S. Federal Law restricts the device to sale by or on the order of a dentist or physician.

Ver. 2506-001

GC

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED by
GC BIOMATERIALS
3737 West 127th Street, Alsip, IL 60803 U.S.A.
1-866-563-7339
www.gcbiomaterials.com

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541
TEL.: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL.: +61 2 9301 8200

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