FOR STERILE COMPONENTS

Prior to use, carefully read the instructions for use and GC Implant Aadva Manual.

GC Implant Aadva® DENTAL IMPLANT SYSTEM

Instruction for use for STERILE COMPONENTS
For use only by a dental professional in the recommended indications.

GC Implant Aadva	Composition
Standard Implant <p>Tapered Implant Short Implant Cover Screw Healing Screw SR Abutment SR Abutment 17° SR Abutment 30° Ball Abutment</p>	<p>Titanium alloy, (Titanium-6%Aluminium-4%Vanadium)</p>

1. INDICATIONS

- Indications for implant treatment:
 - Partial and total edentulism
 - Poor tooth prognosis
 - Improvement of oral-health related quality of life
 - Avoidance of tooth preparation
 - Development of natural tooth appearance
 - Preservation of crestal volume
 - Improvement of denture retention and stabilisation
- Special indications – GC Aadva Short Implant

GC Aadva Short Implant is indicated for use in the mandible with reduced vertical bone dimension that can accommodate GC Aadva Short Implant without impeding the mandibular canal and its content. Short implant is to be placed with the rough portion being submerged, leaving the polished part at the crestal level. In the case of partial bridges, A Ø4.2-mm Short Implant should always be restored with a splinted prosthetic restoration including additional adjacent implant(s). Ø5.2 mm Short Implants can be used for single-unit restorations of cases with good crown-to-implant ratio, low masticatory forces, appropriate intermaxillary distance and absence of parafunctional habits or occlusal diseases.
- Special indications – GC Aadva Narrow Implant

GC Aadva Narrow Implant is indicated for patients with low masticatory forces, appropriate intermaxillary distance and absence of parafunctional habits or occlusal diseases. The following regions present feasible recipient sites:

 - Maxillary lateral incisors
 - Mandibular incisors

If used in the premolar region, it is recommended that it is restored with a splinted prosthetic restoration including additional adjacent implant(s). GC Aadva Narrow Implant is not indicated for use in the posterior and canine region.

2. CONTRAINDICATIONS

- Absolute contraindications:
 - Hypersensitivity to titanium or components of titanium alloy
 - General medical state precluding minor oral surgical procedures (i.e. recent myocardial infarction/CVI, severe blood disorders, bleeding issues, uncontrollable endocrine disorders, recent valvular prosthesis, severe renal disorder, etc.)
 - High dose of radiation and/or chemotherapy
 - Intravenous bisphosphonate therapy
 - High risk for medication-induced osteonecrosis of the jaw
 - Acute inflammation in the recipient region
 - Uncontrolled periodontal disease
 - Mental health status precluding invasive surgery
- Relative contraindications:
 - Low dose radiation therapy in the head and neck region
 - Immunocompromised patient
 - Systemic and metabolic disease and/or treatment impairing healing and bone physiology
 - Oral bisphosphonate therapy
 - Low risk for medication related osteonecrosis of the jaw
 - Incomplete growth and development
 - Inadequate width and height in recipient bone bed or extraction socket
 - Chronic inflammation, benign lesions (i.e. tumors, cysts) or trauma in the recipient region
 - Parafunctional habits (i.e. bruxism, clenching), trismus, unfavourable intermaxillary distance, occlusal diseases
 - Heavy drug, alcohol and/or tobacco use
 - Poor oral hygiene

PREPARATION

- Implant
 - Non-sterile assistant: Open blister pack and drops sterile container on sterile cloth.
 - Steile assistant: Hold container with the cap facing upward and remove the cap.
 - Take out the implant.
- Abutment/Screw

Non-sterile assistant: Peel open non-sterile blister pack and drop the component into sterile titanium bowl.

PRECAUTIONS

- GC Implant Aadva may only be used by dental professionals having received appropriate education and training.
- 1 Follow instructions for use as outlined in the GC Implant Aadva Manual.
- 1 Prior to surgery, a thorough general medical and oral patient examination including radiographic and model analysis should be performed and a treatment plan should be carefully created.
- Components of the GC Implant Aadva system cannot be used in combination with other implant systems.
- 5 Do not use if the sterile container is damaged or previously opened. Do not use the damaged or contaminated product.
- 5 Sterile components are for single-use only. Do not re-sterilize. Reuse or re-sterilization of these components could cause patient infection or failure/loss of osseointegration. Do not use after expiration date.
- Always use sterile instruments specially provided for this system.
- If any sensitivity reactions are experienced, discontinue use of the product and refer to a physician.
- Regular examination and follow-up treatment after surgery is important.

ADVERSE EVENTS

- Failure to or loss of osseointegration.
- Prolonged pain, numbness, changes in sensation.
- Delayed type dental metal allergies.

STORAGE

Store in a cool dry place away from heat, moisture and direct sunlight. Shelf life: as described on the package

Magnetic Resonance Conditional
Non-clinical testing has demonstrated that the Aadva Implant System configuration is MR Conditional. A patient with Aadva Implants and Abutments can be safely scanned in an MR system meeting the following conditions:

- Static Magnetic Field Strength: 1.5 or 3 T
- Maximum spatial field gradient:
 - 1.5 T: 60.37 T/m (= 6037 G/cm)
 - 3 T: 30.19 T/m (= 3019 G/cm)
- Radiofrequency (RF) field exposure:
 - Radiofrequency excitation: Circularly Polarized (CP)
 - Radiofrequency transmit coil: whole-body transmit coil
 - Radiofrequency receive coil type: whole-body receive coil
 - Maximum permitted whole-body averaged specific absorption rate (SAR): Normal Operating Mode, 2 W/kg

Note:
Under the scan conditions defined above, the Aadva Implant system configuration is expected to produce a maximum temperature rise of less than 2.3 °C at 1.5 T and less than 2.6 °C at 3 T after 15 minutes of continuous scanning.

In non-clinical testing, the image artefact caused by one implant extends approximately 12 mm from the device edge on gradient echo image, and 13 mm on spin echo images, at 1.5 T.

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