EN

Prior to use, carefully read the instructions for use.

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GC Fuji IX

GLASS IONOMER FOR ATRAUMATIC RESTORATIVE TREATMENT

GC Fuji IX is a glass ionomer restorative material specially designed for the ART technique (Atraumatic Restorative Treatment). The cavity preparation technique is quite different from normal operative dentistry and the application technique has been modified in order to allow the operator to achieve a high level of success with limited facilities.

ADVANTAGES

- Simple reliable technique for all operators.
- Chemical adhesion to dentin and enamel.
- Continuous fluoride release.
- Adequate placement time and rapid set.
- Excellent biocompatibility.
- Minimal water sensitivity
- High compressive strength
- Proven track record

RECOMMENDED INDICATIONS

- Class I and small Class II restorations in deciduous teeth.
- Class I and small Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.
- Fit and sealing.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.





PROPERTIES ON USAGE

Powder/Liquid ratio (g/g)	3.6/1.0
Mixing time (sec.)	25'-30"
Working time (min., sec)"	2'00"
Net setting time (min., sec)"	2'20"

Test conditions: Temperature (23±1°C), Relative Humidity (50±10%)

DIRECTIONS FOR USE

1) Cavity Preparation
a) Using gloved fingers, isolate the tooth with cotton rolls, thus keeping the treatment area free from saliva.

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
BRIEF REVIEW

Note:
1. Avoid use of GC Fuji IX with patients who have a history of allergy to glass ionomer cements.
2. Do not use GC Fuji IX with patients who show an allergy to the material. In case of allergy, immediately stop using it and advise them to consult a physician.
3. The operator who has a history of allergy to glass ionomer cements should not handle GC Fuji IX. In case of allergy, he should immediately stop using it and consult a physician.
4. Do not allow the mixture to contact the oral tissues or skin. In case of contact, immediately remove the material with absorbent cotton soaked in alcohol and rinse with water.
5. Avoid contact of the mixture with eyes. In case of contact, immediately flush with water and seek medical attention.
6. Clean the nozzle of the liquid bottle occasionally with a damp cotton pad to avoid clogging.
7. If pulp capping is required, a calcium hydroxide such as GC Nu-Cap can be used to protect pulp.
8. After use, tightly close both liquid and powder bottles to prevent exposure to moisture.
9. Do not mix the powder or liquid of GC Fuji IX with any other glass ionomer product.
10. The material is designed for use at temperature of 21-25°C (69.8-77.0°F).
11. This product is for use only by a dental professional.
12. This product should be used only for the applications described in the instructions.

SHADE

Universal shade

PACKAGE


GC Fuji IX 1-1 pack:
1 x 10g powder, 1 x 6g (4.8ml) liquid, 2 x 2g GC COCOA BUTTER, plastic spatula, powder scoop, and a mixing pad (#22)

GC Fuji IX 1-1 mini pack:
1 x 5g powder, 1 x 4g (3.2ml) liquid, 1 x 2g GC COCOA BUTTER, powder scoop

STORAGE

Store in a cool and dark place (4-25°C / 39.2-77.0°F).
Shelf life
GC Fuji IX powder and liquid: 3 years from date of manufacture
GC COCOA BUTTER: 5 years from date of manufacture

Last revised : 12/2013

GC Fuji IX

Glassionmer cho ky thuat tram rintang khong sang chan

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