Prior to use, carefully read the instructions for use.

GC Gold Label Luting & Lining Cement

RADIOPAQUE SELF CURE GLASS IONOMER LUTING & LINING CEMENT

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

- Cementation of metal-based inlays, onlays, crowns, bridges and posts.
- Cementation of high strength (zirconia based) all-ceramic crowns and bridges.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement.

DIRECTIONS FOR USE

Powder / Liquid (g/g)	1,8/1,0
Mixing Time (sec.)	20"
Working Time (23°C / 73°F) (min., sec.)	2'00"
Net Setting Time (37°C / 99°F) (min., sec.)	2'30"
Finishing Time (37°C / 99°F, after seating the restoration) (min., sec.)	4'30"

ISO 9917-1 Glass polyalkenoate cement, Luting

- TOOTH PREPARATION
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared tooth with pumice and water.
 - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).
- RESTORATION PREPARATION
Pre-treat the restoration and handle according to manufacturer's instructions.
- DISPENSING OF POWDER AND LIQUID
 - The standard powder to liquid ratio is 1.8 g/1.0 g (1 level large scoop of powder to 2 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid).
 - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - Hold the liquid bottle vertically and squeeze gently.
 - Close bottles immediately after use.
- MIXING
Dispense powder and liquid onto the mixing pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds. Note:
When mixing smaller amounts, mix rapidly for 15-20 seconds. When mixing larger amounts, divide the powder into two equal portions. Mix the first portion with all of liquid for 5 seconds. Incorporate the remaining portion and mix thoroughly for 15 seconds (total 20 seconds).
- CEMENTATION
 - Coat the internal surface of the restoration with sufficient cement. Seat immediately. At 23°C (73°F) working time is 2 minutes from the start of mixing. Higher temperatures shorten the working time.
 - Maintain moderate pressure on the seated restoration.
 - Start removing excess cement when the excess cement feels rubbery.
 - Finishing can be started 4 minutes 30 seconds after seating the restoration.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39-77°F).

SHADE

Light yellow

PACKAGES

- 1-1 package: 35 g powder, 25 g (20 mL) liquid, double-ended powder scoop, mixing pads (No. 20)
- Bottle of 35 g powder with double-ended powder scoop
- Bottle of 25 g (20 mL) liquid
- Mini pack: 15 g powder, 10 g (8 mL) liquid, double-ended powder scoop

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- This product is not indicated for filling or core build-up.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<https://www.gc.dental/america>
<https://www.gc.dental/sea>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

Last revised: 06/2026

GC

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1708-597-0900
www.gc.dental/america

GC SOUTH AMERICA COMÉRCIO DE
PRODUTOS DENTAIS LTDA
Rua: Maria Cumpati, 441 - Sala 5020 D
Santana, São Paulo, SP 05451L
CEP: 02452-001 CNPJ: 08.279.999/0001-61
TEL: +55-11-2925-0965 RESP. TEC: Enick de Lima CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541
TEL: +65 6546 7588

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