Prior to use, carefully read the instructions for use.

ENGLISH

GC Gold Label Universal Restorative

RADIOPAQUE GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Restoration of primary teeth.
2. Core build-up.
3. Restorations of Class III, V and limited Class I cavities.

CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2.7 / 1.0
Mixing time (sec.)	25"-30"
Working time (min. sec.)	2'00"
Net setting time (min. sec.)	2'20"

Test conditions : Temperature (23+/-1°C), Relative humidity (50+/-10%) ISO 9917-1 : 2003(E) (Glass Polyalkenoate Cement)

1. POWDER AND LIQUID DISPENSING

- a) The standard powder to liquid ratio is 2.7/1.0g. 1 level scoop of powder to 1 drop of liquid.
- b) For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- c) Hold the liquid bottle vertically and squeeze gently.
- d) Close bottles immediately after use.

2. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (total time within 30 seconds).

3. RESTORATIVE TECHNIQUE

- a) Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping use calcium hydroxide.
- b) Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER or GC DENTIN CONDITIONER is recommended to remove the smear layer.
- c) Mix the required amount of cement. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- d) Remove surface moisture, but DO NOT DESICCATE.
- e) Transfer cement to the preparation using a syringe or suitable placement instrument. Avoid air bubbles.
- f) Form the contour, and place a matrix if required during the first 2 minutes of setting.
- g) Setting time is 5 minutes 30 seconds after start of mixing. Once set, remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (Light Cure).

4. FINISHING

- a) One visit treatment :
 - 1) After initial set, correct the contour WITHOUT WATER SPRAY using standard techniques.
 - 2) Apply GC Fuji VARNISH or GC Fuji COAT LC and wait for about 15 minutes.
 - 3) Final finishing and polishing under water spray can now be carried out.
 - 4) Apply a final coat of GC Fuji VARNISH or GC Fuji COAT LC.
- b) Standard method :
24 hours after placement finish and polish the restoration under water spray using standard techniques.

STORAGE

Store at temperature of 4-25°C (39.2-77.0°F). (Shelf life : 3 years from date of manufacture)

SHADE

No. 21 - pale yellow
No. 22 - yellow brown
No. 23 - dark grey
No. 26 - deep brown

PACKAGES

1. 1-1 Package: 15g powder, 10g (8.0mL) liquid, powder scoop, mixing pad (No. 22)
2. Bottle of 15g powder with scoop.
3. Bottle of 10g (8.0mL) liquid.
4. Mini pack*: 10g powder, 7g (5.6mL) liquid, powder scoop

* Mini pack is only available in limited shades (No. 21, 22).

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. DO NOT mix powder or liquid with any other glass ionomer components.

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