Fig.1 / Abb.1

Fig.2 / Abb.2


Fig.3 / Abb.3


Fig.4 / Abb.4


Fig.5 / Abb.5

Fig.6 / Abb.6

Fig.7 / Abb.7

Fig.8 / Abb.8

Fig.9 / Abb.9

Fig.10 / Abb.10

Fig.11 / Abb.11

Fig.12 / Abb.12

Fig.13 / Abb.13

Fig.14 / Abb.14

Fig.15 / Abb.15


Fig.16 / Abb.16

Fig.17 / Abb.17

Fig.18 / Abb.18

Fig.19 / Abb.19

Fig.20 / Abb.20

Fig.21 / Abb.21

Fig.22 / Abb.22


Fig.23 / Abb.23

Fig.24 / Abb.24

**GC Initial™ LiSi Press**

LITHIUM DISILICATE GLASS CERAMICS

For use only by a dental professional in the recommended indications.

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**CONTRA-INDICATIONS**

Not contraindicated for the following patients:

• Severe reduced occlusal dimension

• Reduced cervical and/or labial area by 0.3mm, and the incisal edge by 0.4mm.

**PREPARATION DESIGN**

When designing restorations, the following preparation guidelines and minimum dimensions should be maintained:

Preparation guidelines for all ceramic restorations

• Margins with deep chamfer or rounded shoulder. (fig. 1)

• Avoid having margins in direct occlusal contact with the opposing tooth.

• Antagonist contacts must be taken into consideration.

Occlusal veneer

• Reduce the anatomical shape and observe the stipulated minimum thickness.

• Reduce the incisal crown third – in the occlusal area by approx.1.0mm.

Thin veneer

• If possible, the preparation should be located in the enamel.

• Reduction in the cervical and/or labial area by 0.3mm, and the incisal edge by 0.4mm.

Inlays

• A preparation depth of at least 1.0 mm and an isthmus width of at least 1.0 mm must be observed in the fissure area.

• Do not prepare undercuts.

• Ensure that the walls form an angle of 6° with the long axis of the tooth.

• All internal edges and angles should be rounded.

Onlays

• A preparation depth of at least 1.0 mm and an isthmus width of at least 1.0 mm must be observed in the fissure area.

• Do not prepare undercuts.

• Ensure that the walls form an angle of 6° with the long axis of the tooth.

Crowns in the anterior region

• Occlusal veneer should be at least 1.0mm.

• Reduce the incisal crown third – in the occlusal area by approx.1.5mm.

3-unit bridges

• The maximum acceptable pontic width is different in the anterior and posterior regions.

• The pontic width is determined on the unprepared tooth.

Crowns in the posterior region

• Reduce the anatomical shape and observe the stipulated minimum thickness.

• Width of the shoulder/chamber should be at least 1.0mm.

• Reduce the incisal crown third – in the occlusal area by approx.1.5mm.

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