# *initial*™ LiSi Press

## Technical Manual

# Master all your challenges

Photoshooting: MDT M. Brüsche, Germany  
MDT, S. Maffei, Italy  
MDT Björn Pfeiffer, Germany  
MDT Gábor Filep & DT Max Martini, Germany  
MDT Francesco de Giovanni, Italy

Photoshooting step-by-step: MDT M. Brüsche, Germany

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EN Thank you for choosing GC Initial LiSi Press, lithium disilicate glass ceramics.  
GC Initial LiSi Press is the revolutionary new pressable ceramics that combines unparalleled strength and exceptional aesthetics with the major bonus of a faster, more user-friendly process.  
GC Initial LiSi Press is optimized to be used with the rest of the GC Initial family, including the already proven GC Initial LiSi veneering ceramic, GC Initial IQ SQIN and GC Initial Lustre Pastes NF/ONE.  
This technical manual will give you an impression of how easy it is to get a convincing aesthetic result with minimum workload and excellent physical properties.

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
EN The result is the ultimate combination of strength and aesthetics, making GC Initial LiSi Press perfectly suitable for all types of restorations through all levels of translucency. Critically, HDM technology helps ensure the product remains super stable, without distortion or drop in value, after multiple firings.

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


# 3





## EN Indications For Use DE Indikationen

- Veneers
- Inlays
- Onlays
- Crowns in the anterior and posterior region
- 3-unit bridges in the anterior region
- 3-unit bridges in the premolar region up to the second premolar as the terminal unit
- Crown or splinted crown on top of an implant abutment
- 3-unit bridges up to the second premolar placed on top of an implant abutment

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EN **Contraindications**

GC Initial LiSi Press is not indicated for patients with:

- Severely reduced residual dentition
- Severely reduced occlusal vertical dimension
- Parafunctions (e.g. suffering from bruxism)
- Patients with known allergies to any of the components in GC Initial LiSi Press.

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EN A simplified shade line up with four translucencies to match all your indications, a reduction of inventory & cost and adaptable for a highly esthetic build-up. From physical properties and processing point of view, all type of ingots can be used for all type of indications. Nevertheless, based upon the four translucency groups, following recommendations in the table beneath can guide you to obtain the most esthetic results.

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
EN **High Translucency**

Color line-up that replaces natural tooth enamel. Does not look dark in the mouth. Ideal for small and minimal invasive restorations such as veneers, inlays and onlays. Characterization can be done using GC Initial Spectrum Stains or with GC Initial<sup>TM</sup> Lustre Pastes NF/ONE.

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High Translucency (HT) / Hohe Transluzenz (HT)


EN **Medium Translucency**

V-Shade line-up with warm colors from the GC Initial family of ceramic. Perfect for small restorations such as veneers, inlays and onlays where less transparency and more dentine color is required; or as a full contour restoration in the posterior region. MT restorations are ideal for the staining technique with GC Initial Spectrum Stains or GC Initial<sup>TM</sup> Lustre Pastes NF/ONE. Esthetics can be enhanced with GC Initial IQ SQIN in the microlayering technique or GC Initial LiSi using the cut-back technique.

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Medium Translucency (MT) / Mittlere Transluzenz (MT)


EN **Low Translucency - Stain or Layer**

Low translucency ingots, following the V-shade line-up. Ideal for staining, microlayering with GC Initial IQ SQIN or cut-back layering with GC Initial LiSi.

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
EN **LT-IQ (One Body Concept A, B, C, D or layer)**

Compact color line-up, following the GC Initial IQ - One Body concept. Due to their lower translucency it is perfect to use for the microlayering technique with GC Initial IQ SQIN or the cut-back or layering technique using GC Initial™ LiSi. In the posterior region LT can be used as a full contour restoration and characterized with GC Initial™ Spectrum Stains or GC Initial™ Lustre Pastes NF/ONE.

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
EN **Medium Opacity**

Thanks to strong fluorescence, a life-like sense of color can be reproduced when veneering with GC Initial™ LiSi. Ideal for layering or the cut-back technique. MO ingots are ideal to mask discolored teeth or titanium abutments. In the posterior region, MO can even be used as a full contour solution, and characterized with GC Initial™ Spectrum Stains or GC Initial™ Lustre Pastes NF/ONE.

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Medium Translucency (MT)

Low Translucency (LT)

Low Translucency - One Body Concept (LT-IQ)

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
\*Shades can be achieved by applying GC Initial<sup>TM</sup> Spectrum Stains, GC Initial<sup>TM</sup> Lustre Pastes NF/ONE or veneering.

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# 5

## EN Preparation Design DE Präparationsdesign

EN When designing all ceramic restorations, the following preparation guidelines and minimum dimensions should be maintained.


Basic guidelines for all ceramic restorations:

- No angles or edges
- Prepare margins with deep chamfer or rounded shoulder
- The indicated dimensions should match the minimum thicknesses of the GC Initial LiSi Press restorations
- The margin line of the preparation should not be in contact with the antagonist
- Antagonist contacts must be taken into consideration

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


### Occlusal veneer Okklusales Veneer

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



### Onlay


- A preparation depth of at least 1,0 mm and an isthmus width of at least 1,0 mm must be observed in the fissure area.
  - Do not prepare undercuts.
  - Ensure that the cavity walls form an angle of 6,0° with the long axis of the tooth.
  - All internal edges and angles should be rounded.
  - Provide at least 1,0 mm of reduction in the cusp areas.
- DE
- Eine Präparationstiefe von mindestens 1,0 mm und eine Isthmus-Breite von mindestens 1,0 mm muss im Fissurenbereich eingehalten werden.
  - Präparieren Sie keine Unterschnitte.
  - Stellen Sie sicher, dass die Kavitätswände einen Winkel von 6,0° mit der Längsachse des Zahns bilden.
  - Alle inneren Kanten und Winkel sollten abgerundet werden.
  - In den Höckerbereichen muss eine Reduzierung von mindestens 1,0 mm vorgesehen werden.

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

### 3-unit bridges 3-teilige Brücken

- EN
- The maximum acceptable pontic width is different in the anterior and posterior region. The pontic width is determined on the unprepared tooth.
  - In the anterior region (up to the canine), the pontic width should not exceed 11 mm.
  - In the premolar region (up to the second premolar), the pontic width should not exceed 9 mm.

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
## EN Framework Design DE Gerüstdesign

**EN** The design of the lithium disilicate glass ceramic framework contributes to the longevity and durability of all-ceramic restorations. A well designed framework provides a high quality final result, clinical success and patient satisfaction. The framework design instructions below for the staining, cutback and layering technique should be followed in order to obtain satisfying esthetical. Furthermore, should they be taken in consideration to meet the preparation guidelines.

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



## EN 6.1 Framework design for staining technique

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
## EN 6.2 Framework design for cutback and layering technique

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### EN Model preparation

Fabricate the working model with removable segments as usual using GC Fujirock EP. Coat dies with a die hardener. Use a sealer and a spacer up to 1 mm from the cervical margin.

#### Modelling /contouring

Use an organic wax or a resin that burns out without leaving any residues.

Respect minimum wall thicknesses (refer to page 17 - 19) and design framework according to processing technique (staining , cut-back or layering technique).

Exact contouring, particularly in the preparation area is recommended to avoid time consuming and risky fitting procedures.

For full contour restorations it is recommended to consider occlusal relief, as the glazing process slightly modifies the vertical dimensions.

#### Sprueing

Use a silicone ring system suitable for lithium disilicate press ceramic technique (Ø13 mm)

The wax restorations are sprued onto the ring base former using wax.

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


**NOTE**

- Always attach the sprues in the direction of the flow of the ceramic material and at the thickest part of the wax pattern so that smooth flowing of the viscous ceramics during pressing is achieved.
- If only one object is invested and pressed in a furnace, a second short (blind) sprue must be placed.
- Make sure to keep the space noted in the figure for the wax pattern to assure stability during the pressing procedure.
- Weigh the wax restorations including the sprues correctly. The definitive wax weight is 0.4g (min.) – 0.8g (max.) for 1 ingot.

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**TIP**

- When multiple investment rings are preheated together, place a mark onto the ring base with wax indicating the shade.
- Margins should be facing outwards.

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## EN Investing DE Einbetten

**EN** The investing is carried out with GC LiSi PressVest, a specially developed carbon-free phosphate bonded speed investment for press-ceramic techniques offering you an easier, faster and safer processing. Hydrofluoric acid treatment and alumina blasting are no longer required for removing the reaction layer. Please refer to the IFU of GC LiSi PressVest.

**GC LiSi PressVest offers you a number of advantages:**

- Easy removing of the reaction layer
- Easy devesting
- High fluidity and long working time
- Precise expansion control
- Optimized for quick heating technique
- Adapted to multiple press ceramic techniques
- Reduction of operating time
- Reproducible fit due to no use of hydrofluoric acid
- Reduction of operating time
- Accurate & detailed reproduction
- Optimal fit

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- EN SR Liquid contains a high concentration of surface refining agent for easier removing of the reaction layer.  
Spray just once (1x) on the internal / external surface of the wax pattern.  
Keep a distance of 15cm for an optimal dispersion of the spray.  
Remove immediately residual SR Liquid in / on the wax object with air blowing.


- NOTE**
- Shake the spray bottle well.
  - Residue and excess SR Liquid may cause rough surface. Be sure that all mist is dried out before investing.
  - Do NOT use SR Liquid for multi-sprue technique such as bridges. It may emphasize the fused ceramic part.
  - Wipe off the excess SR Liquid from the former before use to avoid contamination.

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- EN Apply a slight layer cocoa butter on the surface of the ring base and ring gauge.  
This will assure an easy removal after setting.

Carefully place the silicone ring onto the ring base without damaging the wax objects.

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- EN Expansion can be adjusted by diluting the liquid by diluting the mixing liquid with distilled water.

Basic rule: the higher the liquid concentration, the higher the expansion.

The concentration can be freely adapted based on the working experience of the technician.

Please refer to the IFU of GC LiSi PressVest.

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- EN Exact powder/liquid measurement is necessary to obtain stable results. Use adapted measuring equipment such as electric balance, liquid measuring cup, pipette.

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EN Investing  
DE Einbetten

EN 8.4 Mixing

Pre-mix powder and liquid thoroughly for 30 seconds by hand with a spatula.  
Be sure all powder is wetted out by the liquid to give a uniform mixture.  
Always pour liquid first in the mixing bowl, before adding the investment powder.

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EN Mix for 60 seconds under vacuum (320-420rpm).  
Always use clean mixing bowl and check vacuum level.  
Insufficient vacuum leads to differences in fit and air-bubbles of the pressed objects.

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EN 8.5 Working time

GC LiSi PressVest has super high flow during 6 minutes after start of mixing at room temperature (23°C/73°F).  
Working time depends on the temperature of powder & liquid and on room temperature.  
Higher temperatures shorten the working time.

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## EN 8.6 Investing

Carefully fill the cylinder with investment material just below the edge.

GC LiSi PressVest is very fluid, so strong vibration is not necessary and not advisable.

From the moment when the ring is totally filled, stop vibration immediately and do NOT touch the investment until set.

**TIP** Apply a small quantity of investment material onto the gauge. Then carefully place the gauge onto the silicone ring with a rotating movement.

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
### EN 8.7 Setting time

Leave to set for 20 minutes from the start of mixing.

After setting, remove the mould base former and ring base former. Carefully push the investment ring out of the elastic cylinder. Smooth the bottom side by using a sharp instrument. Check the 90° angle and stable position.

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EN A longer bench time is possible up to a maximum of 180 minutes from the start of mixing.

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- EN Place the investment ring with the funnel facing down in a pre-heated burn-out furnace at 850°C/1560°F and keep it for minimum 45 minutes

**NOTE** This material is exclusive to quick heating. Please do NOT apply conventional step-heating technique. Quick heating up at the maximum possible rate results in a relatively constant expansion of the investment, avoiding the formation of cracks and preserving the strength of the investment material.

In case of placing several investment rings in the furnace at one time, extend the holding time depending on the number of investments to prevent surface defect. Due to aggressive burning out, do not open the furnace during heating-up.

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### EN TIPS

- The ring should be placed in the furnace between 20 minutes to 180 minutes from the start of mixing. Longer bench setting could lead to higher risk of cracks.
- Preheat the burn-out furnace to 900°C, once the ring is placed in the furnace, lower the temperature to 850°C.
- Do not mix metal- and press ceramic rings in the same furnace. The metal rings excessively lower the temperature of the furnace. At maximum, four 100 g rings or two 200 g rings can be placed at the same time in the furnace.
- The holding time should be a least 45 minutes at 850°C.
- Avoid lowering the temperature by opening the ring furnace during the holding time for burnout.

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### EN 10.1 Pressing procedure

1. Remove the investment ring from the preheating furnace immediately after completion of the preheating process.
2. Insert a GC Initial LiSi Press ingot into the investment ring with printed side facing up.
3. Place a plunger into the investment ring. Now place the investment ring in the center of the press furnace table.
4. Start the dedicated press program.

#### NOTE

It is recommended to use one-way press plungers. Do not preheat the plunger or ingot. Only one ingot per investment ring may be used for pressing. Please follow the ceramic press furnace manufactures' instructions for use. The ingot should be inserted as quickly as possible to prevent the ring from cooling down, before placing it into the press furnace. Avoid rapid cooling after the end of press cycle.

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
Dekema Austromat® is a registered trademark of DEKEMA Dental-Keramiköfen GmbH, Freilassing, Germany

Programat® are registered trademarks of Ivoclar Vivadent AG.

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## EN 10.1 Pressing programs

### Important:

The above mentioned pressing parameters are only guidelines and therefore always need to be adjusted to the pressing furnace and its correct functionality.

These pressing parameters can only be used as guidelines.

Most important is to obtain the right pressing result. We refer to Presstemperature calibration chart.

**For other furnaces than those mentioned below, please consult your furnace manufacturer or dedicated GC salesforce.**

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
EN This chart visually helps you to evaluate your press result and adapt your press temperature accordingly. Raise or lower press temperature gradually with 5°C and evaluate. To obtain a perfect press result, it is necessary that also the press cylinder has the correct temperature.

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# 11

## EN Devesting DE Ausbetten

- EN**
1. After cooling, mark the length of plunger on the investment ring.
  2. Cut along the marking with a suitable cut-off disk.
  3. Carefully break into the 2 sections.
  4. Rough devesting is carried out by blasting with glass beads at 4 bar pressure.
  5. Fine devesting is carried out by blasting with glass beads at 2 bar pressure.
  6. Carefully cut off and finish the sprues with a suitable cut-off disk and smoothen the area.

**NOTE**

- Do not use aluminum oxide for sandblasting. Using GC LiSi PressVest investment guarantees the reduction of the reaction layer and can be accordingly completely removed by blasting with glass beads.
- You do NOT need to put pressed objects in hydrofluoric acid.
- Avoid overheating during the cutting off of the sprues.

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
# 13

- EN Glazing, staining & layering
- DE Glasieren, Bemalen und Schichten

- EN For glazing, staining and layering, use our dedicated GC Initial Lustre Pastes NF/ONE, GC Initial IQ SQIN, GC Initial LiSi veneering ceramic & GC Initial Spectrum Stains for the best esthetic results.

Sandblast the surface of the restoration under low pressure (50 microns Alu-Oxid, 1.5 bar) before applying the GC Initial™ Lustre Pastes NF/ONE, GC Initial IQ SQIN or GC Initial™ LiSi.

Please refer to the IFU of Initial IQ Lustre Pastes NF/ONE, GC Initial IQ SQIN and GC Initial LiSi for glazing, staining, and layering.

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EN 13.1 GC Initial IQ Lustre Pastes NF/ONE

DE 13.1 GC Initial IQ Lustre Pastes NF/ONE

EN The Initial IQ Lustre Pastes NF/ONE 3-dimensional ceramic stains bring colour deepness and lifelike translucency in all your ceramic restorations and are in a ready-to-use consistency. The use of these pastes is very versatile: due to their optimal match to both "low" and "high" CTE ceramics, they can be used for the complete GC Initial ceramic range.

The Lustre Pastes NF/ONE are based on fine ceramic particles and can be applied in a thicker layer where they will exhibit excellent vitality and a very natural glaze to your restorations.

On lithium disilicate restorations, the GC Initial IQ Lustre Pastes NF/ONE can be used for


- the glazing and characterization of monolithic, full anatomical, lithium disilicate restorations
- the glazing and characterization of the GC Initial LiSi veneering material
- the wash firing of the lithium disilicate structure, providing an immediate internal characterization, before microlayering with GC Initial IQ SQIN or layering with GC Initial LiSi ceramics.

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EN 13.2 GC Initial LiSi  
DE 13.2 GC Initial LiSi

- EN GC Initial LiSi is a feldspar veneering ceramic adapted to the light-dynamics of lithium disilicate glass-ceramic frameworks guaranteeing a highly aesthetic and durable restorations for long-term patient satisfaction. GC Initial LiSi is characterized by a precisely adapted coefficient of thermal expansion (CTE), low firing temperatures and maximum stability even after multiple firings. It will impress you with its excellent aesthetics on lithium disilicate glass-ceramic frameworks, particularly in terms of their optical dynamics.

This veneering material can be used both in the individual layering and in the cut-back technique following the proven Initial ceramic concept - One shade system - One layering technique. If you want to go even further, full anatomical frameworks can be enhanced by staining using the GC Initial Lustre Pastes NF/ONE.

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EN Crowns which have been pressed and finished according to manufacturer specifications: tooth 16 monolithic, tooth 15 partially to be veneered, tooth 14 completely to be veneered.

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EN Important note: The surfaces to be glazed with GC Initial IQ Lustre Pastes must first be sandblasted using  $\text{Al}_2\text{O}_3$  at 1 bar and subsequently cleaned with steam cleaner.

First coat the sandblasted surfaces of the "not to be veneered" areas extensively with Lustre Paste L-N.

To achieve the desired V-shade colour, rub the body Lustre A to D more or less intensively into the neutral lustre.

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EN Coat all fully anatomical surfaces with the L-N lustre and the body Lustre A-D. For the Wash Firing sequence of the portions to be veneered, individualization can be done by applying the body Lustre A-D or the Effect Lustre in a thin layer covering the entire surface.

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
EN The glaze-fired crowns with a wide variety of GC Initial IQ Lustre Paste applications can be individualized with GC Initial LiSi veneering ceramic.

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EN More information about step-by-step veneering can be found in the GC Initial LiSi Technical manual

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EN More information about the use of GC Initial Spectrum Stains be found in the Technical Manual.

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### EN 13.5 GC Initial<sup>TM</sup> FIRING FOAM

### DE 13.5 GC Initial<sup>TM</sup> FIRING FOAM

EN When firing, do not heat or cool the restorations quickly. Rapid change in temperature could break the material. When firing, proper furnace tray (honeycomb tray) and support pin in combination with the GC Initial FIRING FOAM should be used. The GC Initial FIRING FOAM ensures a stable fixing of objects onto the ceramic firing tray.

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## EN Cementation DE Zementierung

### EN Clear and simple luting protocol

Lithium disilicate glass ceramics restorations can be adhesively or conventionally luted, depending on the indication and the practitioner's preferences.

The GC Luting Guide will assist you in selecting the best alternatives for each indication.

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
### EN 14.1 G-CEM ONE

### DE 14.1 G-CEM ONE

EN G-CEM ONE is a truly universal self-adhesive resin cement for all indications and all type of restorations. Depending on the preparation, you could choose from either a simplified self-adhesive cementation or full adhesive application.

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- EN 1. Apply hydrofluoric acid gel (5-9%) for 20 seconds to the inner surfaces of the restoration.  
2. Rinse with water spray or an ultrasonic cleaner and dry.  
3. Condition the etched surfaces with a silane coupling agent such as Ceramic Primer II or G-Multi PRIMER and allow it to dry.

\* Phosphoric acid (35-37%) can be also used for the purpose of cleaning the surface, preferably scrubbing with a microbrush for 10-15 seconds.

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EN Cementation  
DE Zementierung

EN In case the preparation is non-retentive, an adhesive resin cement (such as G-CEM ONE with dedicated bonding) is preferred. Download our GC Luting Guide for step-by-step instructions for each cementation option.

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EN Based upon a dedicated shape via an STL-file you can mill / print yourself shade tabs that can be pressed in GC Initial LiSi Press. In this way you easily your own LiSi Press shade tabs. These 3D printed, castable shade tabs are ready to invest & press fitting into the existing GC Initial Shade guides.

Special thanks to MDT Stefan Roozen (Austria) for the idea and images of this individual LiSi Press shade tabs.

Contact your local GC representative for the STL of the dedicated shape to make your own castable shade tabs or to request a sample of the castable shade tab.

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Case by MDT. S. Roozen, Austria / Fall von MDT. S. Roozen, Österreich

---

**EN Medical Device Regulation Requirement**

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en) as well as to our internal vigilance system: [vigilance@gc.dental](mailto:vigilance@gc.dental)  
In this way you will contribute to improve the safety of this product.

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RESPONSIBLE MANUFACTURER IN CANADA  
GC AMERICA INC.  
3737 W. 127th Street, Alsip, IL 60803 U.S.A.

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
Tel: +1-708-597-0900  
<https://www.gc.dental/america/>

GC Germany GmbH  
Seifgrundstraße 2, D-61348 Bad Homburg  
Tel. +49.61.72.99.59.60, Fax. +49.61.72.99.59.66.6  
info.germany@gc.dental, <https://www.gc.dental/europe/de-DE>

GC ITALIA S.r.l.  
Via Calabria 1, I-20098 San Giuliano Milanese  
Tel. +39.02.98.28.20.68, Fax. +39.02.98.28.21.00  
info.italy@gc.dental, <https://www.gc.dental/europe/it-IT>

GC UNITED KINGDOM Ltd.  
Coopers Court, Coopers Court, Newport Pagnell, UK-Bucks. MK16 8JS  
Tel. +44.1908.218.999, Fax. +44.1908.218.900  
info.uk@gc.dental, <https://www.gc.dental/europe/en-GB>

GC FRANCE s.a.s.  
8 rue Benjamin Franklin, F-94370 Sucy en Brie Cedex  
Tel. +33.1.49.80.37.91, Fax. +33.1.45.76.32.68  
info.france@gc.dental, <https://www.gc.dental/europe/fr-FR>

GC IBÉRICA - Dental Products, S.L.  
Edificio Codesa 2  
Playa de las Américas, 2, 1º, Of. 4, ES-28290 Las Rozas, Madrid  
Tel. +34.916.364.340, Fax. +34.916.364.341  
comercial.spain@gc.dental, <https://www.gc.dental/europe/es-ES>

GC AUSTRIA GmbH  
Tallak 124, A-8103 Gratwein-Strassengel  
Tel. +43.3124.54020, Fax. +43.3124.54020.40  
info.austria@gc.dental, <https://www.gc.dental/europe/de-AT>

GC Europe NV - Benelux Sales Department  
Researchpark Haasrode-Leuven 1240  
Interleuvenlaan 33, B-3001 Leuven  
Tel. +32.16 74.18.60  
info.benelux@gc.dental, <https://www.gc.dental/europe/fr-BE>

GC EUROPE N.V. - East European Office  
Siget 19b, HR-10020 Zagreb  
Tel. +385.1.46.78.474, Fax. +385.1.46.78.473  
info.eeo@gc.dental, <https://www.gc.dental/europe>

GC NORDIC AB - Finnish Branch  
Lemminkäisenkatu 46, FIN-20520 Turku  
Tel. +358.40.900.07.57  
info.finland@gc.dental, <https://www.gc.dental/europe/fi-FI>

GC NORDIC AB - GC Nordic Danish Branch  
Scandinavian Trade Building  
Gydevang 34-41, DK-3450 Allerød  
Tel. +45 51 15 03 82  
info.denmark@gc.dental, <https://www.gc.dental/europe/da-DK>

GC NORDIC AB  
c/o Lundin Revisionbyrå  
Erik Dahlbergsgatan 11B, SE-411 26 Göteborg  
Tel. +46.768.54.43.50  
info.nordic@gc.dental, <https://www.gc.dental/europe/sv-SE>

GC AUSTRIA GmbH - Swiss Office  
Zürichstrasse 31, CH-6004 Luzern  
Tel. +41.41.520.01.78, Fax. +41.41.520.01.77  
info.switzerland@gc.dental, <https://www.gc.dental/europe/de-CH>

GC AUSTRALASIA DENTAL PTY LTD  
1753 Botany Rd, Banksmeadow NSW 2019, Australia  
Tel: +61-2-9301 8200, Fax: +61-2-9316 4196

GC SOUTH AMERICA  
Rua Heliadora, 399, Santana - São Paulo, SP, BRASIL  
CEP: 02022-051 - TEL: +55-11-2925-0965 - CNPJ: 08.279.999/0001-61  
RESP. TÉC: Erick de Lima - CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.  
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541  
Tel: +65 6546 7588