Prior to use, carefully read the instructions for use.

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FR

SOLARE X

LIGHT-CURED COMPOSITE RESTORATIVE

For use only by a dental professional in the indications for use.

INDICATIONS FOR USE

1. Direct restorative for Class I, II, III, IV, V cavities.
2. Direct restorative for wedge-shaped defects and root surface cavities.
3. Direct restorative for veneers and diastema closure.

CONTRAINDICATIONS

1. Direct pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

PRODUCT DESCRIPTION

SOLARE X is a light-cure, radio-opaque restorative material to be used intra-orally and classified as a Type 1 and Class 2 (Group 1) per ISO standard 4049. This material has a radiopacity equivalent to 2.0-2.5 mm of aluminium (dentin = 1 mm, enamel = 2 mm). The particle size of inorganic fillers range is 0.02-2.9 µm. The total amount of inorganic filler is approximately 65 vol%.

COMPOSITION

Fluoro-alumino-silicate glass, dimethacrylate, feldspar, contrast agent, multifunctional monomer, silicon dioxide, initiator, pigment, stabilizer


DIRECTIONS FOR USE

1. Shade Selection
Clean the tooth with pumice and water. Shade selection should be made prior to isolation.
2. Cavity Preparation
Prepare cavity using standard techniques.
Note:
For pulp capping, use calcium hydroxide.
3. Bonding Treatment
For bonding SOLARE X to enamel and / or dentin, use a light-cured bonding system such as SOLARE UNIVERSAL BOND, or G-BOND or GC Fuji BOND LC. Follow manufacturer's instructions.
4. Placement of SOLARE X
Dispensing from a syringe
Remove syringe cap and dispense material onto a mixing pad. Place the material into the cavity using a suitable placement instrument. After dispensing, screw syringe plunger anticlockwise by a half to full turn to release residual pressure inside the syringe. Replace cap immediately after use.
Note:
1. Material can be applied in a single layer to achieve aesthetic restorations using Standard shades. For details, refer to the Clinical Hints.
2. Material may be hard to extrude immediately after removing from cold storage. Prior to use, leave to stand for a few minutes at normal room temperature.
3. After dispensing, avoid too long exposure to ambient light. Ambient light can shorten the manipulation time.

Clinical Hints

- a) In the case of small cavities
Restore using a one shade technique. In most cases the use of one Standard shade alone will be sufficient. In cases where a higher degree of translucency is needed, one of the Outside special shades can be selected.
- b) In the case of large and/or deep cavities
In most cases a multi shade layering technique will give the best aesthetic results. To block out shine through from the oral cavity or to mask discolored dentin, place X-AO2 and continue to build up with a Standard shade. For optimal aesthetics use an Outside special shade as the final composite layer. In the case of deep posterior cavities, a flowable composite such as SOLARE Flo or a glass ionomer cement such as GC Fuji LINING LC (Paste Pak) or GC Fuji IX GP or Gold Label 9 EXTRA can be used on the cavity floor instead of an Opaque shade. See also Examples of Clinical Applications.

Examples of Clinical Applications (Clinical Hint a, b)

5. Contouring before Light Curing
Contour using standard techniques.
6. Light Curing
Light cure SOLARE X using a light curing unit. Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Irradiation Time and Effective Depth of Cure

Shade	Irradiation time	High power LED (more than 1200 mW/cm²)	10 sec.	20 sec.
X-A1, X-WT	X-A2, X-A3	Halogen / LED (700 mW/cm²)	20 sec.	40 sec.
X-A1, X-WT			3.0 mm	3.5 mm
X-A2, X-A3			2.5 mm	3.0 mm
X-A3.5, X-AO2			2.0 mm	3.0 mm

The effective wavelength range of each dental curing unit must cover 450-480 nm.

Note:

1. Material should be placed and light cured in layers. For maximum layer thickness, please refer to the table.
 2. Lower light intensity may cause insufficient curing or discoloration of the material.
7. Finishing and Polishing
Finish and polish using diamond burs, polishing points and discs. To obtain a high gloss, polishing pastes can be used.

SHADES

Standard shade: X-A1, X-A2, X-A3, X-A3.5

Inside special shade: X-AO2

Outside special shade: X-WT

A, AO shades are based on Vita® * Shade.

*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C / 39.2-77.0°F) away from high temperatures or direct sunlight.

PACKAGES

1 syringe of 3.7 g (2 mL)

Accessories (available separately)

- a. Shade guide
- b. Mixing pad (No.14B)

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Take care to avoid ingestion of the material.
4. Wear plastic or rubber gloves during operation to avoid direct contact with air inhibited resin layers and prevent possible sensitivity.
5. Use a protective light shield or similar protective eye wear during light curing.
6. When polishing the polymerised material, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
7. Do not mix with other similar products.
8. Avoid getting material on clothing.
9. In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
10. Do not use SOLARE X in combination with eugenol containing materials as eugenol may hinder SOLARE X from setting.
11. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
12. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

<http://www.gceurope.com>

They can also be obtained from your supplier.

CLEANING AND DISINFECTING:

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at Regulatory.gce@gc.dental

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

Last revised: 12/2023

GC

MANUFACTURED by
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 100

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541
TEL: +65 6546 7588

GC AUSTRALASIA DENTAL PTY. LTD.
1753 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

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