Prior to use, carefully read the instructions for use.

EN

TEMPSMART™

DUAL-CURED COMPOSITE FOR TEMPORARY CROWN AND BRIDGE

For use only by a dental professional for the following indications.

RECOMMENDED INDICATIONS

- Fabrication of temporary crowns, bridges, inlays, onlays and veneers
- Fabrication of long-term temporary restorations

CONTRAINDICATIONS

Avoid use of this product in patients with known methacrylate monomer and methacrylate polymer allergies.

DIRECTIONS FOR USE

Fabricating temporary restorations

- Taking Impression
 - Take impression using alginate or silicone impression material. Alternatively, temporary can be made using a laboratory-prepared thermoforming matrix.
- Block out preparation undercuts.
 - Clean tooth preparation with water spray and lightly dry.
 - Lightly lubricate with separating medium, such as petroleum jelly on prepared teeth or resin core.
- Dispensing

NOTE:

Before dispensing first time from a new syringe or after a long interval between uses, bleed pastes from the syringe.

 - Remove cap of syringe by rotating 1/4 turn counter-clockwise. Attach a Mixing Tip by aligning V shaped notch on rim of the mixing tip with V shaped notch between syringe barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of mixing tip 1/4 turn clockwise.
 - Discard first centimeter of extruded material before use.
 - Extrude material directly into impression.

- NOTE:
- Used mixing tip can be left on syringe to serve as a storage cap until next use. For each application, place a new mixing tip.
- Fabricating a temporary restoration

Fabricate a temporary crown and bridge material according to timing shown below. Time indicated starts after mixing paste.

Option 1: Fabricating in the mouth

0'00"~0'45"	Place TEMPSMART into impression and reposition in mouth.
0'45"~2'00"	Setting in mouth.
2'00"~2'30"	Remove impression from teeth.
	If necessary, tack cure each tooth unit with curing light for 1~2 seconds.

Final cure

	Remove temporary restoration from mouth and final cure each tooth unit with light curing unit (Plasma Arc (2000 mW/cm), G-Light (1200 mW/cm), Halogen / LED (700 mW/cm)) for 5 seconds. Irradiate occlusal surface. If chemical cure only, keep temporary restoration in 23°C/74°F until 5'00".
--	---

After final cure

	Completely remove oxygen inhibition layer from surface of temporary restoration with alcohol.
--	---

Option 2: Fabricating on the model

0'00"~0'45"	Place TEMPSMART into impression and reposition on model.
0'45"~3'00"	Setting on model.

Final cure

	Remove impression from the model and final cure with an extra oral light curing oven for 1 minute.
	Remove temporary restoration from model.
	If chemical cure only, keep temporary restoration in 23°C/74°F until 5'00".

After final cure

	Completely remove oxygen inhibition layer from surface of temporary restoration with alcohol.
--	---

Caution!

Due to high polymerization temperature created, do not irradiate temporary restoration for more than 2 seconds. Do not final cure on tooth.

- NOTE:
- Above processing times are applicable for products stored and processed at a temperature of 23° C / 74° F and at 50% relative humidity. Higher temperatures shorten these times, while lower temperatures extend them. The size of the restoration also has impact on setting times: as size of restoration increases, setting takes longer.
 - When irradiating temporary crown and bridge, irradiate occlusal surface.
 - When final curing temporary crown and bridge with curing light, irradiate one tooth from both ends first to prevent deformation.
 - Finishing
 - Remove excess material using fine carbide bur.
 - Check occlusion and if necessary, correct shape by using a resin material such as TEMPSMART, G-ænial Flo, G-ænial Universal Flo or UNIFAST Trad, or grind to obtain proper occlusion.
 - Polish surface of temporary restoration using silicone polisher and goat's hair wheel. Alternatively, apply a coating agent such as OPTIGLAZE or characterize with OPTIGLAZE Color following manufacturer's directions.

Caution!

Do not breathe in polishing dust. Use personal protective devices such as facemask and safety glasses.

NOTE:

 - When using a resin material such as TEMPSMART, G-ænial Flo, G-ænial Universal Flo or UNIFAST Trad for custom-shaping and shape correction, usually there is no need for a bonding agent. Completely remove saliva contaminated part using a carbide bur. If necessary, apply a ceramic primer such as CERAMIC PRIMER II to surface of temporary restoration.
 - When using TEMPSMART for custom-shaping and shape correction, light cure mixing paste with curing light for 10 seconds.
 - Cementing

Cement using temporary cement and follow manufacturer's instructions

NOTE:

Do not use eugenol cement as it will prevent polymerization of resin.

Repairing / Modifying temporary restorations

 - Abrade surface of temporary restorations using carbide bur.
 - If necessary, apply a ceramic primer such as CERAMIC PRIMER II and dry with air.
 - Place a resin material such as TEMPSMART, G-ænial Flo, G-ænial Universal Flo or UNIFAST Trad and cure according to manufacturer's instructions.

NOTE:

Completely remove saliva contaminated part using a carbide bur.

SHADES

A1, A2, A3, A3.5, B1 and BW

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25° C / 39.2-77.0° F) away from high temperatures or direct sunlight.

PACKAGES

Refills:

- Cartridge 48 mL (60 g) (1), MIXING TIP II SSS (16)
- Cartridge 10 mL (12.5 g) (1), GC Automix Tip Regular (10)

Mixing tip package:

 - MIXING TIP II SSS (60)
 - GC Automix Tip Regular (20)
- Dispenser package
 - CARTRIDGE DISPENSER II (1)
 - GRADIA CORE DISPENSER GUN (1)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate operating field from skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Take care to avoid ingesting material.
- Do not mix with components of other products.
- Mixing tips are for single use only. To prevent cross contamination between patients, do not reuse the mixing tip. The mixing tip cannot be sterilized in an autoclave or chemiclave.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn. Wearing gloves during use avoids direct contact with air inhibited resin in order to prevent possible sensitivity.
- Wear protective eye glasses during light curing.
- When polishing polymerized material, use dust collector and wear a dust mask to avoid inhalation of dust.
- Avoid getting material on clothing.
- In case of contact with unintended areas of tooth or prosthetic appliances, remove with instrument, sponge or cotton pellet before light curing.
- Do not use TEMPSMART in combination with eugenol containing materials as eugenol may delay or prevent TEMPSMART from setting.
- Do not wipe mixing tip with excessive alcohol as this may cause contamination of the paste.
- Dispose of all wastes according to local regulations.
- In rare cases product may cause sensitivity in some people. If any such reactions are experienced, discontinue use of the product and refer to a physician.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered health care grade infection control product according to regional / national guidelines.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:

http://www.gcamerica.com

They can also be obtained from your supplier.

Last Revised: 06/2026

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torinamiya-cho, Kasugai, Aichi 486-0844, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

RESPONSIBLE MANUFACTURER IN CANADA:
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900
www.gc.dental/america

GC SOUTH AMERICA COMÉRCIO DE PRODUTOS DENTAIS LTDA
Rua: Maria Curupaiti, 441 - Sala 5020 D
Santana, São Paulo, SP, BRASIL
CEP: 02452-001 CNPJ: 08.279.999/0001-61
TEL: +55-11-2925-0965 RESP. TEC: Erick de Lima CRO/SP 100.866

GC ASIA DENTAL PTE. LTD.
5 Tampines Central 1, #06-01 Tampines Plaza 2, Singapore 529541

GC AUSTRALASIA DENTAL PTY. LTD.
1763 Botany Rd, Banksmeadow, NSW 2019, Australia
TEL: +61 2 9301 8200

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