

COE-SOFT™ PRINT

Resilient Denture Liner

The advanced, cost-effective solution for denture relining



Since 1921

100 years of Quality in Dental



Combining trusted quality with modern innovation

COE-SOFT™ PRINT is specially designed for ease, comfort, and precision. Whether for temporary liners in upper or lower dentures or post-surgical care, this resilient, self-curing material meets diverse clinical needs with efficiency.

COE-SOFT™ PRINT simplifies the process, enhances post-surgical care, and delivers outstanding results. With its advanced features and thoughtful design, it is the go-to solution for modern denture relining needs.

Features & Benefits:

- **3D Material Compatibility:** Bonds effortlessly to 3D printed materials for enhanced adaptability and precision
- **Cost-Effective:** Offers premium results at a competitive price, making it an economical choice for both practices and patients
- **Perfect for Post-Surgical Applications:** Offers gentle support and soft comfort, making it ideal for sensitive cases following surgical procedures
- **Quick Mixing Time:** Ready to use in just 30 seconds, ensuring fast preparation and application
- **Mild Mint Flavor:** Provides a light, refreshing flavor for added patient satisfaction
- **Efficient Polymerization:** Cures in or out of the mouth in approximately 15 minutes for streamlined workflows
- **Temporary Durability:** Ensures reliable performance, lasting up to three months for temporary relines
- **Soft Comfort:** Maintains the trusted soft-touch experience patients appreciate.



10041721 COE-SOFT™ PRINT Professional Package

Contains: one powder (156g), one liquid (163mL; 169g), one bonding agent (12g, 12.7mL), one measuring cup, one mixing cup; one vial, one spatula, one lubricant.

10041722 COE-SOFT™ PRINT Liquid Refill

Contains: One refill bottle (163mL; 169g).

10041723 COE-SOFT™ PRINT Powder Refill

Contains: One refill bottle (156g).

10041724 COE-SOFT™ PRINT Bonding Agent

Contains: One refill bottle (12.7mL; 12g).