




G-aenial™ Universal Flo

UNIVERSAL INJECTABLE COMPOSITE RESTORATIVE

THE FIRST TRULY UNIVERSAL INJECTABLE COMPOSITE



Results of G-aenial Universal Flo Evaluation
at 1 Year from The Dental Advisor



TOUCH of GENIUS™

- Superior polish and high gloss retention¹
- Lack of post-operative sensitivity
- Extraordinary wear resistance
- Strength through filler innovation - for Classes I, II, III, IV and V
- Wears like a traditional composite, yet handles like a flowable



The First Truly Universal Injectable Composite

3M

Truly Universal - Indicated for:



Class I



Class II



Class III



Class IV




Class V

True innovation in filler size and composition, in filler silanation and in manufacturing techniques, have given rise to a remarkable composite resin. G-ænial Universal Flo is the first truly universal Injectable Composite, a material that clinicians can place with ease for beautifully adapted, beautifully finished, fast, aesthetic and tough composite restorations.


Restored using shades A03, A4, A3 and AE.




Dr. S. Koide, Nagano, Japan



G-ænial Universal Flo - Features 200-nm strontium glass fillers homogeneously dispersed for high flexural strength and wear resistance



Filtek Supreme Ultra Flow[†] - (Nano composite). Typically contains spherical nano clusters made up of agglomerated nano fillers of varying sizes



Clearfil Majesty[†] - (Micro Hybrid). Typically contains a range of glass and nano fillers of varying sizes

400222 **G-ænial[™] Bond & G-ænial[™] Universal Flo Kit**
Contains: One G-ænial Bond Kit (5 mL), 1 Syringe A2 G-ænial Universal Flo and 20 Dispensing Tips.

G-ænial[™] Universal Flo Refills

Contains: One Syringe and 20 Dispensing Tips.
004202 A1 004203 A2 004204 A3 004205 A3.5 004206 A4 004216 AE
008468 AO1 004211 AO2 004212 AO3 004207 B1 004208 B2 004209 B3
004213 BW 004210 C3 004214 CV 004215 JE

G-ænial[™] Flo Refills

Contains: One Syringe and 20 Dispensing Tips.
004194 A1 004195 A2 004196 A3 004197 A3.5 004198 A4 004199 AO2
004200 AO3 004201 CV

004635 **Metal Dispensing Tips III (x30)**

004634 **Plastic Dispensing Tips III (x30)**

GC

¹Data on file.
[†]Not a registered trademark of GC America Inc.