



# EXPERIENCE EFFICIENT LIGHT TECHNOLOGY

## The Light & The Light 405

Cordless LED Curing Light

### Two Lights - One Great Result Every Cure!

- Concentrated light provides deep and efficient polymerization
- Reduce eye fatigue by removing unnecessary dispersed light
- Easy-to-read OLED screen
- Variety of curing modes and time settings available
- Lightweight handpiece with 330° head twist
- Interchangeable, long lasting battery
  - Over 400 cures in one charge (high mode)
  - Over 1,000 cures in one charge (low mode)




**The Light** operates from 440 nm - 490 nm (peaking at 460 nm) to meet all standard light-curing needs. Its light intensity per mode is: High (up to 1500 mW/cm<sup>2</sup>), Low (up to 650 mW/cm<sup>2</sup>), Soft Start (up to 1500 mW/cm<sup>2</sup>) and Pulse (up to 1500 mW/cm<sup>2</sup>).


**The Light 405** operates from 405 nm - 490 nm to meet all standard light curing needs, with the same intensity modes as The Light. The Light 405 includes a dual wavelength processor designed to meet all standard operator light curing needs, including that of OPTIGLAZE™ color, GC's Light-Cured Characterization Coating for Composite and Acrylic Indirect Restorations.



The Light (uniform concentrated light increases polymerization performance)



Competitor LED (insufficient light dispersion)



## Synergize with OPTIGLAZE™ color

The Light 405 has a dual wavelength used to cure all light cured restorative materials, as well as OPTIGLAZE™ color, our light-cured characterization coating for composite and acrylic indirect restorations.

With OPTIGLAZE™ color, you can now change shades, shift value, characterize, create effects and glaze all your resin type prosthetics in no time at all. Best of all, offer your dentists and patients a more beautiful and esthetic restoration that will last due to GC's innovative new homogeneously dispersed nano-filler technology.



### 402400 The Light (Operator)

Contains: One hand piece, one charger, one battery, one user manual and one box of 100 disposable covers.

### 402405 The Light (Operator & Laboratory)

Contains: One hand piece, one charger, one battery, one user manual and one box of 100 disposable covers.

### 402402 Replaceable Battery

Contains: One battery.

### 402403 Disposable Covers

Contains: 100 replaceable disposable covers.

### 008408 OPTIGLAZE™ color Set

Contains: All color shades (2.6 mL each), disposable palette, flat and round brush and brush holder.

Scan Here For a Product Demo

