


Speed Meets Performance:



CEREC Workflows with CERASMART® 270

Same-day dentistry, even faster.



CEREC Software 5.3.4 introduces the new Super Fast Grinding process — delivering significantly faster restorations without changing the workflow you already trust.

Paired with CERASMART® 270, this update delivers faster, more efficient chairside workflows — without firing and without additional equipment.

CERASMART 270 is ideally suited to this enhanced workflow.

As a force absorbing, hybrid ceramic material, it mills efficiently, finishes easily, and supports high-quality, aesthetic restorations — making it a smart choice for fast, same-day dentistry.

Indication	Super-Fast	Fast	Fine
Crown	~4 min	7:14 min	11:54 min
Inlay	~1:40 min	6:17 min	10:18 min



Since 1921
100 years of Quality in Dental



Less waiting. More efficiency.

Trusted CERASMART® 270 results —
now delivered faster.

What this means for your practice:

- Faster grinding times with CEREC 5.3.4
- Quality: Smaller sprue for composites (2mm → 1.4mm), new “purple” diamond tools for hybrid/composite, high restoration quality.
- No firing required with CERASMART® 270
- Reduced chair time and improved practice efficiency
- Consistent, reliable results using a proven and CEREC-validated material



New “purple” diamond tools for hybrid/composite only. Used for Super-Fast, Fast and Fine mode.

[VIEW BROCHURE](#)

[SHOP CERASMART 270](#)

Upgrade your CEREC software and fast-track restorations with CERASMART® 270.



Since 1921
100 years of Quality in Dental