# Essentia<sup>®</sup> Universal shade Ultimate simplification for all posterior restorations





# Essentia<sup>®</sup> Universal shade THE solution for all posterior restorations

Essentia® Universal shade offers complete flexibility and simplicity for all posterior restorations.

Just one shade, with amazing blendability, lets you create aesthetic restorations with ease.

Essentia<sup>®</sup> Universal shade is now available in three full-strength viscosities, so you can select the perfect solution for each clinical indication. Choose between:



# One shade. Absolute simplicity.



When GC introduced Essentia<sup>®</sup>, it opened the door to a world of simpler restorative dentistry. Now that door is thrown wide open with Essentia<sup>®</sup> Universal shade.

## Three full strength viscosities

#### **Essentia**®





A direct approach was selected for this challenging post-endodontic case



Restoration with everX Posterior™ as a dentin core and Essentia Universal as the outer shell



Six-month recall showing good shade integration with the surrounding dental tissues

#### Essentia<sup>®</sup> LoFlo





Air abrasion of the non-carious lesions



Restoration with Essentia® LoFlo Universal



Post-operative view after polishing

#### Essentia® HiFlo





Preparation of the small Class I cavities



Restoration using Essentia® HiFlo Universal



Dr Bojidar Kafelov, Bulgaria

Based on a mix of

Based on homogeneously

dispersed 200nm strontium glass fillers

**Essentia**<sup>®</sup> is a full-strength paste composite with a simplified one shade concept. Its Universal shade displays a great blending effect together with a packable-consistency **paste viscosity**, for easy placement and shaping.

pre-polymerised and ultra-fine fillers

**Essentia**<sup>®</sup> **LoFlo** is a full-strength injectable composite with **low flowability**, optimal when a thixotropic material is needed – for instance for cervical restorations. It promises minimal wear, ensuring long-lasting restorations, even on the occlusal surface.

**Essentia® HiFlo** is a full-strength injectable composite with **high wettability**, perfect for easy adaptation to preparations and flowing into narrow cavities. It offers great handling properties, together with superior strength and a high radiopacity of over 300% (1mm %Al), facilitating the diagnosis of secondary caries.



# Essentia<sup>®</sup> Universal shade

## Essentia<sup>®</sup> Universal shade refill

Essentia® Syringe, U (1 x 2mL) Essentia® LoFlo Syringe, U (1 x 2mL, 20 dispensing tips)

Essentia® ge, HiFlo Syringe, 20 U (1 x 2mL, 20 tips) dispensing tips) Essentia® Unitips, U (15 x 0.16mL)



## Essentia<sup>®</sup> Universal shade kit



## **G-Premio BOND**

G-Premio BOND unit doses, 50 pieces (0.1mL) G-Premio BOND Bottle Refill 5mL G-Premio BOND 3-bottle Pack 5mL liquid (3)



essent



a **, 'C.C', '** 



GC Australasia Dental Pty Ltd 1753 Botany Road Banksmeadow NSW 2019 Australia T: +61 2 9301 8200 E: info@gcaustralasia.com www.gcaustralasia.com