

Intuitive Digital Workflows for Your Daily Challenges



The Intraoral
Scanner from GC
- Aadva™ IOS 100 P
based on AI Technology



GC

Aadva IOS 100 P - addressing your daily practice needs

IOS 100 Portable



Portable

IOS 100 P is a compact system - only a handpiece, cradle and laptop

Anti-fog heating

Scans are easily maintained with the anti-fogging function in the cradle and the powder and spray-free technology

Light handpiece

IOS 100 P has impressively small and light handpiece

Easy File Transfer

USB, local network and Aadva Xchange

No Calibration Needed

Trouble-free and reliable results every time you scan

Open System

Open STL and PLY files offering full compatibility with third party software and hardware

Visual and Audio Feedback

Gain full control of your scanning procedure


through its modular, intuitive, case-management solutions

General scanning features

- AI Technology filters the soft and moving tissues to guarantee a smooth scanning
- Scan, pause, check and rescan anytime you want
- Automatic bite alignment
- Simple and intuitive user interface



Ortho Module



- Customised workflow for your ortho case follow up
- Dedicate PDF file to assist treatment planning
- Possibility to record both central relation (CR) centric occlusion (CO) bites
- Integrated ortho model maker, hollowed out (American-cut) for optimal 3D printing
- Integrated third party communication possibilities (VDDS)
- Scan files compatible with various aligner systems
- Patented impression scan workflow

Prostho Module



- Customised workflow for your crown and bridge work
- Automatic margin line for your preparations
- Single implant scan
- Simple order sheet with necessary information for labs
- Automated import of files to exocad



Aadva™ Xchange

All-in-one approach for digital case-management. The cloud-based platform supports communication with the GC production and service teams as well as with registered external dental labs.



Trainings

We educate and empower dentists and dental lab technicians with innovative solutions in digital dentistry.



Support

GC Europe's customer support for its digital dental technology is at industry-leading level.



LinkedIn



YouTube



GC EUROPE N.V.

Head Office
Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info.gce@gc.dental
www.gceurope.com

GC NORDIC AB

c/o Lundin Revisionbyrå
Erik Dahlbergsgatan 11B
SE-411 26 Göteborg
Tel. 0768 54 43 50
Info.nordic@gc.dental
<http://www.gceurope.com>