

<b>VITA® Classic</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A4</b>	<b>B1</b>	<b>B2</b>	<b>B3</b>	<b>B4</b>	<b>C1</b>	<b>C2</b>	<b>C3</b>	<b>C4</b>	<b>D2</b>	<b>D3</b>	<b>D4</b>
Opaak	AO1	AO2	AO3	AO3	AO3	AO1	AO2	AO3	AO3	AO2	AO3	AO3	AO3	AO1	AO2	AO3
Põhitoon	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A4</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A2</b>	<b>A3</b>	<b>A3</b>	<b>A3.5</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>
Email	JE	JE	AE	AE	AE	JE	JE	AE	AE	JE	AE	AE	AE	JE	JE	AE

Mitmehilise tehnika vastavus

G-ænial® **A'CHORD**

Email



Opaak



Põhitoon



Email



Põhitoon



Email



**Süstal**  
(1 x 2.1ml/4.0g)  
**Kapslid**  
(15 x 0.16ml/0.33g)

10006779

10006753



**Põhitoonid** A1  
**VITA® Classic** A1 B1 D2

10006780

10006754



A2  
A2 B2 C1 D3

10006781

10006755



A3  
A3 B3 C2 C3 D4

10006782

10006756



A3.5  
A3.5 B4 C4

10006783

10006757



A4  
A4

## G-ænia<sup>®</sup> A'CHORD

### Opaagid

BOW



10006790

10006764

AO1



10006787

10006761

AO2



10006788

10006762

AO3



10006789

10006763

### Standard toonid

BW



10006786

10006760

A5



10006784

10006758

A6



10006785

10006759

### Emailid

JE



10006791

10006765

AE



10006792

10006766

**Süstal**  
(1 x 2.1ml/4.0g)  
**Kapslid**  
(15 x 0.16ml/0.33g)