

Spectrum Stains



SPS-1
Ivory White



SPS-2
Melon Yellow



SPS-3
Signal Yellow



SPS-4
Light Terracotta



SPS-5
Savannah Gold



SPS-6
Caramel



SPS-7
Camel



SPS-8
Olive



SPS-9
Brown



SPS-10
Dark Walnut



SPS-11
Coral



SPS-12
Denim Blue



SPS-13
Twilight



SPS-14
Deep Orange



SPS-15
Sunset



SPS-16
Midnight



SPS-17
Dark Brown



SPS-18
Illusion 1



SPS-19
Illusion 2




SPS-20
Illusion 3

SPS-17, 18, 19, 20 are not part of the Spectrum Stain Set (art.nr. 877094)



This chart is intended as a reference guide only.


Step by step




Coat the entire surface with Spectrum Glaze



SPS-15
Application of sunset



SPS-14
Application of deep orange



SPS-13
Application of twilight



SPS-16
Application of midnight



SPS-1
Application of ivory white



Before (L): Ceramic crown with individual micro-texture
After (R): Perfectly stained and glazed crown with complete preservation of the micro-structure.



SPS GL



SPS GL-FLU

Glaze Fluo is not part of the Spectrum Stain set

