

As easy as 1-2-3

Repair Kit from GC

Smart products
and guidelines
for ALL your intra-oral
repair procedures

GC

The GC Repair Kit, chairside repair made easy...

A **complete set** of all products and guidelines needed to cover all your **routine intra-oral repairs** in an **easy and smart way**.



Metal-based indirect restorations (precious and non-precious)



Indirect restorations made of **glass-ceramics, zirconia, alumina or hybrid ceramics** (e.g. CERASMART)



Composite restorations (direct or indirect)

What's in it for you? Three components to fit all repairs!

Ceramic Primer II

Following the zero compromise approach, we suggest the use of a **separate silane primer** (Ceramic Primer II) together with our adhesive G-Premio BOND, for a durable adhesion to **glass ceramics** and **hybrid ceramics**.



GRADIA PLUS Opaque

This **light-cured opaquer** efficiently masks dark backgrounds, to ensure high aesthetics in case of porcelain fused to metal restorations.

G-Premio BOND

A **universal bonding** which offers **optimal adhesion to all substrates** thanks to 3 different monomers.

4-MET

for a reliable adhesion to **enamel & dentin**

MDP

for an excellent adhesion to **enamel & dentin, zirconia, alumina and non-precious metals**

MDTP

to bond to **precious metals***



For the step-by-step, consult also the GC Restorative Dentistry Guides!

Article Number	Description
901139	GC Repair Kit: Ceramic Primer II (5 unitdoses), G-Premio BOND (5 unitdoses), GRADIA PLUS OAI (1 syringe 2.0mL)



*MDTP: Methacryloyloxydecyl dihydrogen thiophosphate to bond to precious metals

GC EUROPE N.V.

Head Office
Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info@gceurope.com
http://www.gceurope.com



GC UNITED KINGDOM Ltd.

16-23, Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info@uk.gceurope.com
http://uk.gceurope.com