



Cytrans[®] Elashield

Information for patients

1.) What is Cytrans Elashield made of?

Cytrans Elashield is a bioresorbable membrane made of (100%) L-lactide- ϵ -caprolactone copolymer.

2.) What is expected lifetime of the device?

Cytrans Elashield life-time duration is 6 months.

3.) Are there any necessary follow-ups required after initial procedure?

As with any surgical procedure, follow-up visits should be scheduled in a manner to avoid post-operative infection of the site.

Date of issue: 18/03/2026

Version: 1