



# Cytrans<sup>®</sup> Granules

## Information for patients

### 1.) What are Cytrans Granules made of?

Cytrans Granules are a synthetic bone graft material composed of 100% carbonate apatite, a substance with a composition similar to the inorganic mineral component of human bone.

### 2.) What is the expected lifetime of the device?

Cytrans Granules are naturally resorbed over a period of up to 24 months.

### 3.) Are there any necessary follow-up visits required after initial procedure?

As with any surgical procedure, follow-up visits should be scheduled to avoid post-operative infection of the site.

*Date of issue: 02/03/2026*

*Version: 1*