

Vario Healing Screw

GC Tech. Europe

Improves Soft Tissue Management




Since 1921
100 years of Quality in Dental



Vario: for Improved Soft Tissue Management

Aadva, from GC, offers a comprehensive line of healing screws for soft tissue management. Regardless of the clinical situation, Aadva healing screws provide improved aesthetics. They are available in various anatomical shapes allowing for the formation of an emergence profile adapted to prosthetic restorations. Vario healing screws are 100 % compatible with all Aadva implants.



100 % Compatible with Aadva



Colour Coding

Aadva healing screws are colour-coded based on implant diameter. Identification is quick and simple, reducing the potential of mishandling.

We give a guarantee to our clients that we deliver practical solutions based on GC standards

Aesthetic Abutments

The aesthetic abutment line is the latest addition to GC's Aadva portfolio of cemented prosthetic solutions. With the reduction of its emergence profile height to 0,4 mm EPH, aesthetic outcome and integration of the prosthetic abutment are improved, even in cases where the soft tissue level is critical. The newly added healing screw Vario Low, benefits from a treatment process that is both user-friendly and effective, thanks to its perfectly matched EPH to the new aesthetic abutment line.

Ready Abutments



Smart Abutments



Part Numbers

Healing screw

Type	Designation	P/N
● Narrow	Healing Screw Narrow Vario Low	810295
	Healing Screw Narrow Vario Medium	810296
	Healing Screw Narrow Vario High	810297
● Regular	Healing Screw Regular Vario Low	810298
	Healing Screw Regular Vario Medium	810299
	Healing Screw Regular Vario High	810300
● Wide	Healing Screw Wide Vario Low	810301
	Healing Screw Wide Vario Medium	810302
	Healing Screw Wide Vario High	810303

Ready Abutment

● Narrow	Ready Abutment Narrow Hexed Est. 4.5	810304
● Regular	Ready Abutment Regular/Wide Hexed Est. 4.5	810305
● Wide		

Smart Abutment

● Narrow	Smart Abutment Narrow Est. D 4.5	810306
● Regular	Smart Abutment Regular/Wide Est. D 4.5	810307
● Wide		
● Narrow	Smart Abutment 15° Narrow Est. D 4.5	810308
● Regular	Smart Abutment 15° Regular/Wide Est. D 4.5	810309
● Wide		



For more information
please contact:

GC Tech.Europe GmbH

Harkortstr. 2
D-58339 Breckerfeld
Germany
Tel.: +49 2338 801980
Fax: +49 2338 801985
info@gctech.eu
<https://www.gc.dental/gctech/en>