

Torque recommendations



▀ Abutment screw, mono-block abutment and ball abutment: 20 Ncm



▀ SR abutment prosthetic screw: 10 Ncm



▀ Locator abutment: 30 Ncm

Abutment Type	Torque recommendations
Ball Abutment	20 Ncm
Locator™ Abutments	30 Ncm
Prep Abutment	20 Ncm
Provi Abutment	20 Ncm
Smart Abutment	20 Ncm
SR Abutment	20 Ncm
SR Abutment single-unit	20 Ncm
Ti Base Abutment	20 Ncm
UCLA Abutment	20 Ncm
Universal Hybrid Abutment	20 Ncm
SR abutment prosthetic screw	10 Ncm

*Locator™ is a registered trademark of Zest Anchors, LLC