














# デュラデント 前歯

## 上顎



形態番号	COMBINATION	全横径	中切歯横径	中切歯長径
C3		41.0	7.5	9.8
C4		44.0	8.0	10.4
C5		46.0	8.5	11.0
C6		48.5	9.0	11.7

形態番号	COMBINATION-LONG	全横径	中切歯横径	中切歯長径
CL3		41.5	7.5	10.4
CL4		44.5	8.0	11.2
CL5		46.5	8.5	11.8
CL6		49.5	9.0	12.6

形態番号	▲ 尖形 TAPERING	全横径	中切歯横径	中切歯長径
T3		41.0	7.5	9.4
T4		43.5	8.0	10.0
T5		46.0	8.5	10.5

形態番号	● 卵円形 OVOID	全横径	中切歯横径	中切歯長径
O3		40.5	7.5	9.5
O4		43.0	8.0	10.0
O5		45.5	8.5	10.7

## 上顎

形態番号	■ 方形 SQUARE	全横径	中切歯横径	中切歯長径
S3		41.0	7.5	9.5
S4		43.5	8.0	10.0
S5		46.0	8.5	10.7

## 下顎

### STANDARD TYPE

形態番号		全横径	中切歯横径	中切歯長径
3		31.5	4.7	9.2
4		33.5	5.0	9.9
5		35.5	5.3	10.4
6		37.0	5.5	10.8

### LONG TYPE

形態番号		全横径	中切歯横径	中切歯長径
L3		31.5	4.7	10.0
L4		33.5	5.0	10.7
L5		35.5	5.3	11.3

中切歯長径(mm)はカラーを除く。単位は全てmmです。