

**Reliable,
Easy,
Efficient**


No rush, full patient
compliance




Optimax
Fix

G2

Next Generation Developed for Reliable Class II Correction



Optimax Fix G2 is the next generation of Optimax Fix. It is easy-to-use and available in 4 pre-adjusted sizes. Featuring ball and socket hinges for improved comfort, it requires minimum fixation without sacrificing performance.



Due to its design, Optimax Fix G2 offers several advantages:



Easy-to-use

Thanks to 4 pre-adjusted sizes and measurements made directly in the mouth, the system is quickly positioned.




Easy to wear

Reduced size, Optimax Fix G2 minimises patient discomfort and optimum compliance



Comfortable

Maxillary and mandibular ball and socket hinges allow full functional movement for the patient providing optimal comfort



Efficient

Optimax Fix G2 requires minimum fixation:

- Simply insert into the upper molar root tube distally and cinch
- The mandibular arm of the appliance is engaged into the lip bumper tube for stability of the mandibular arch whilst the anterior is engaged and free sliding.

The Optimax Fix G2 is only active when the mouth is closed. No mandibular incisors rotation, No molar intrusion, Maximum control !