



Catalog N°7

# IPR, Lab, Miscellaneous



**Experts in Orthodontics**  
GC Orthodontics Europe GmbH  
[www.gc.dental/ortho](http://www.gc.dental/ortho)



# IPR, Lab, Miscellaneous

Interproximal reduction Dentasonic .....	213
Multi FA .....	217
Laboratory .....	218
Smile Rite Thermoforming .....	218
Extraoral Traction .....	218
Miscellaneous .....	219
Cheek Retractor .....	219
Cheek Retractor with vaccuming system NOLA .....	219



## Interproximal reduction

For safe, precise surface reduction with minimal enamel damage, our professional materials for MI-IPR, reduce inter-proximal spaces with minimal deterioration.



### Profi Kit 2

Profi Kit 2 Box includes single and double-sided DiaStrips, for inter-proximal space reduction, proximal finishing and polishing in orthodontics, aligners and restorative dentistry. Profi Kit 2 also features a new water-cooled oscillating handpiece with a 4:1 ratio and a second head with an alternative movement.

Description	P/N
4 DiaStrips, 2 sides • 4 DiaStrips 1 side • 2 Polishing PoliStrip • 2 polishing burs • 1 water colled hand piece, 4:1 • 1 tête oscillante • 1 measuring gauge set: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm • 1 extracting tip • 1 Spray	55-BD02-0000



### Basic set nwc (no water cooled)

Basic reduction set of 4:1.

Description	P/N
Reduction: 4:1 (not cooled) • 1 assortment of 4 stripes • 1 measuring gauge set: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm • 1 extracting tip	55-BD04-0000



## Basic Set Aligner

Basic set aligner includes single and double-sided DiaStrips, for inter-proximal space reduction, proximal finishing and polishing in aligners and restorative dentistry. Basic Set Aligner also features a new water-cooled oscillating handpiece with a 4:1 ratio with an alternative movement.

Description	P/N
1 DiaSaw 0.1mm, 0.2 mm • 1 DiaStrip, 0.3 mm, 0.4 mm • 3 PoliStrips 15 µm 2 sides • 1 Twist water-cooling HP, 4:1, STROKE 1,1mm • 1 Measuring gauge with 6 thicknesses: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5 • 1 extraction tip • 1 Spray nozzle • Worldwide Warranty.	55-BD06-0000

## Basic Set Ortho

Basic Set ortho Kit includes single and double-sided DiaStrips, for inter-proximal space reduction, proximal finishing and polishing in aligners and restorative dentistry. Basic Set Aligner also features a new water-cooled oscillating handpiece with a 4:1 ratio with an alternative movement.

Description	P/N
1 DiaSaw 0.1mm , 0.2 mm • 1 DiaStrip, 0.2 mm, 0.3 mm, 0.4 mm • 1 DiaStrips polyester 0.5 mm 2 sides • 2 PoliStrips 15 µm • 1 Twist water-cooling HP, 4:1, STROKE 1,1mm • 1 reciprocating Head • 1 Measuring gauge with 6 thicknesses: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5 • 1 extraction tip • 2 years warranty.	55-BD07-0000



## Double-sided DiaStrip

Double-sided DiaStrip oscillating diamond files, allows for inter-proximal reduction on both sides at the same time. Used for proximal finishes and orthodontic polishing.

Description				
Double-Sided	DiaStrip	DiaStrip	DiaStrip	DiaStrip
P/N 5/pk	55-FDM0-5215	55-FDM0-5225	55-FDM0-5240	55-FDM0-5260
P/N 12/pk	55-FDM0-1215	55-FDM0-1225	55-FDM0-1240	55-FDM0-1260
Thickness	0.1mm	0.15mm	0.2mm	0.3mm

## DiaSaw

Diamond oscillating DiaSaw files. Opens contact points safely and quickly (without soft tissue damage) and reduces tooth size in orthodontics.



Double-Sided	DiaSaw	DiaSaw
P/N 5/pk	55-SDM0-5315	55-SDM0-5325
Thickness = +/- 0.03mm	0.1mm	0.15mm

## DiaStrip single side

Single side	DiaStrip left	DiaStrip right
Assortment of 5 pcs	55-FDML-5100	55-FDMR-5100

## Set



Double-Sided	DiaStrip Set 5	DiaSaw Set 5	Aligner Set
P/N	55-FDM0-5200	55-SDM0-5300	55-FDA0-5300
Grain µm:	1x1,0.15,0.2,0.3mm +1PoliStrip	1x1,0.15,0.2mm, +2 PoliStrip	1x1,0.15,0.2mm, +2 PoliStrip
Coating side:	2 sides	2 sides	2 sides
Refill pack:	5 pieces	5 pieces	5 pieces

## PoliStrip

Made from flexible resin, these strips are for ultimate proximal polishing on both sides in just 1 passage. PoliStrips are sold in blister packs of 12 (2x6) with 6 0.10 mm strips and 6 0.15 mm strips. 55-FDRA-1200



### DS-Debonding Polisher

Easily and safely remove orthodontic adhesive residue following bracket debonding. Due to the hardness and abrasiveness of materials used with a rotary bur or a metal instrument, there is always a risk when debonding of damaging the natural tooth structure. Our DS-Debonding Polisher features micro particles deposited on a synthetic rubber matrix that safely removes adhesive residue in minutes.

Description	P/N
Set of 2X3 polisher cones and cup	55-D000-6B00
Set of 6 polisher cones	55-D000-6200
Set of 6 polisher cup	55-D000-6100

### Measuring Gauge

Thickness: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5 mm  
55-G000-1000



### Inox SteryStrip

Stainless steel sterilization set for 6 polishers.  
55-TDM0-6000



### Multi FA

Multi Functional appliance for CL II and CL III. One appliance that fit 90% of cases

Description	P/N
Multi FA Regular	92-5001-0000
Multi FA Short	92-5002-0000
Multi FA Long	92-5003-0000
Multi FA for buccal tubes	92-5004-0000
Multi FA crimpable stops	92-5005-0000





## Laboratory



### Smile Rite Thermoforming

All Smile Rite Thermoforming plastic sheets, regardless of thickness, are pre-dried and packaged in aluminium bags of 25 pieces. The water absorption by the material has been eliminated as well as the bubble problem during thermoforming has been resolved.

Hard	Round 120 mm	Round 125 mm	Square 125mmX125mm
1,5 mm (x 50)	—	92-1015-R125	—
1,0 mm (x100)	92-1010-R120	92-1010-R125	92-1010-S125
0,80 mm (x100)	92-1008-R120	92-1008-R125	—

Soft	Round 120 mm	Round 125 mm	Square 125mmX125mm
1,0 mm (X100)	—	92-1110-R125	92-1110-S125

### Extraoral Traction

These traction devices are made from washable hypoallergenic materials. Comfortable, in vivid colors, they are easily accepted by young patients.



#### Neck Strap with Safety Release

Description	P/N
Blue/Pink 450 g	55-4450-0029
Blue/Yellow 450 g	55-4450-0028
Black 600 g	55-4600-0027

#### Head Cap with Safety Release



Description	P/N
Head Cap Denim 450 g	55-5450-0030
Head Cap Color 600 g	55-5600-0035
Head Cap Denim 600 g	55-5600-0030

## Miscellaneous



### Cheek Retractor

Description	P/N
Opac, Automatic Cheek Retractor, Full View With Extension	55-0240-0012

### Cheek Retractor with vaccuming system NOLA



Description	P/N
Large, white	55-4100-0000
Narrow, red	55-4101-0000
Tongue pusher large by 2	55-4102-0000
Tongue pusher narrow by 2	55-4103-0000
Rétracteur large, white	55-4104-0000
Pipes with connector "Y" white by 4	55-4105-0000
Connector «Y» by 4	55-4106-0000
Silicon pipes par 10	55-4107-0000

