

Based on aerospace applications, the ivory-colored aesthetic-coated Albaflex was subject to intense testing to guarantee its superior behaviour during usage. Albaflex's surface finish limits friction and retention of dental plaque.



Characteristics

- Permanent ivory-colored coating
- Coating does not influence mechanical properties



Aesthetic Archwire **Albaflex**

For a more discreet orthodontic treatment



Need convincing?

Watch the abrasion test here.



Experts in Orthodontics



GC Orthodontics Europe GmbH

www.gcorthodontics.eu

Albaflex Aesthetic Archwires

Albaflex Ivory-colored NiTi Archwires



Albaflex's ivory coating provides an extremely discreet NiTi archwire solution while maintaining its mechanical properties.

		
	References	
Section	Maxillary	Mandibular
.017"	77-C0UC-0017	77-C0LC-0017
.020"	77-C0UC-0020	77-C0LC-0020
.017"x.017"	77-C0UC-1717	77-C0LC-1717
.020"x.020"	77-C0UC-2020	77-C0LC-2020



Albaflex Stainless Steel Ivory-colored Archwires

Albaflex's ivory-colored coating enables perfect discretion for stainless steel archwires while maintaining their mechanical properties.


		
	References	
Section	Maxillary	Mandibular
.017"x.017"	78-C0UC-1717	78-C0LC-1717
.020"x.020"	78-C0UC-1725	78-C0LC-1725
.017"x.025"	78-C0UC-2020	78-C0LC-2020



Albaflex Ivory Spacer

The Albaflex spacer with its outstanding dimensional stability, porosity-free surface and ivory color will convince you of its ability to maintain spaces, protect mucous membrane and stiffen the archwire or make it even more discreet from its very first usage.

Description	Universal
Spacer (ivory-colored) .036", lg 180 mm, 3 pcs	55-0200-0040



99-77EA-1300 REV01-09/2019

