


Access

Standard metal bracket




Access, manufactured by Tomy® in Japan, is an affordable bracket created using MIM (Metal Injection Molded) technology. All information (3 dimensions) are integrated into the bracket. The system features torque in its anatomically laser-milled base, designed to provide retention as reliable as that of traditional mesh.



Features

- Made of injection-molded metal
- Torque in base
- Color ID
- Long axis line
- Anatomical base for each tooth
- In/Out compatible between medium and mini



We recommend using GC Ortho Connect™ for bonding your brackets. This is a one step system with the primer included in the composite.



Combination

Legend mini/medium brackets can be used in conjunction with Chic. brackets to a complete aesthetic result.



Access is manufactured by TOMY Inc. in Yumoto, Japan. Tomy Brackets have been used by orthodontists worldwide since the introduction of their first bracket system in December 1961.

Today, Tomy is recognized as a reputable manufacturer of Orthodontic brackets. Available in MBT Type and Roth Type.

Bracket & Archwire Combination

Combining Access with our Initalloy, BIO-ACTIVE, and BIO-EDGE thermal NiTi archwires. Rhodium-coated archwires are also available for a discreet orthodontic bracket system.

