




GC Ortho Connect™

Light Cure Orthodontic Adhesive

One step system,
with the primer included
in the composite



GC Ortho Connect™ 90-1100-0000 (2x 2,9g)



GC Ortho Connect does not require an additional primer, as it is already contained in the adhesive. This simplifies the bonding process and reduces costs.

Its viscosity ensures reliable and precise fastener placement.

The ergonomic syringe allows precise dosing. Its cap helps to avoid the glue from drying out, minimizing material loss and providing a practical and economical solution.


Adhesive residues can be detectable under black light, aiding you more effective cleaning.

Benefits


- Primer included in the composite
- Fluorescent
- Viscosity
- Resealable applicator tip for precise bonding

High resistance

GC Ortho Connect offers versatility for optimal bonding of metal brackets and tubes, ceramic brackets and plastic brackets.



source: GC Europe NV



source: GC Europe NV



Recommended etching

GC Ortho Etching Gel 90-1400-0000
2x 4,8g + 40 applicator tips

- Good viscosity for precise application
- Remains in place until rinsed off

