



## Advantages

- Integrated primer
- Fluorescence
- Ideal viscosity
- Resealable and precis applicator tip

# GC Ortho Connect™

Light Cure Orthodontic  
Adhesive One Step System




**Experts in Orthodontics**

GC Orthodontics Europe GmbH  
[www.gcorthodontics.eu](http://www.gcorthodontics.eu)


# GC Ortho Connect™ 90-1100-0000 (2 syringes of 2,9 g each)



When using GC Ortho Connect, no primer is necessary, since it is already incorporated into the adhesive, thus creating a simplified bonding process that reduces costs.



The adhesive's ideal viscosity insures reliable and precise bracket placement.




An ergonomic syringe guarantees accurate dosage, providing a practical and economic adhesive solution. The syringe's cap prevents dry out and thus decreasing potential waste.




Excess adhesive is detected by a black light, which aids in fast and efficient cleanup.

## Highly Resistant



### Excellent Bond Strength

Featuring superior flexural strength - fundamental in supporting mastication pressure.



### Reduced Bond Failure

The low elasticity component of GC Ortho Connect, compared to other composites in the industry, allows for high impact resistance without bond failure.

GC Ortho Connect provides reliable and versatile bonding on ceramic and metal brackets and buccal tubes.



## Recommended Etching

GC Ortho **Etching Gel** 90-1400-0000  
2 syringes (4,8 g each) 40 dispensing tips

- Ideal viscosity for precise application
- Stays in place until rinsing

